



SCINTILLATION LIGHT DETECTION IN LARGE LIQUID ARGON TPCS


Ben Jones,
University of Texas at Arlington



Large liquid argon TPCs are central to the US neutrino physics program




DUNE



MicroBooNE Event Zero:

(half design field)

Run 1147 Event 0. August 6th 2015 16:59



MicroBooNE Event Zero:


(half design field)

Run 1147 Event 0. August 6th 2015 16:59

μ BooNE

200 cm

200 cm






MicroBooNE Event Zero:

(half design field)


Run 1147 Event 0. August 6th 2015 16:59

μ BooNE






MicroBooNE simulations





Visibility map




Simulated μ BooNE performance

Plots from
FERMILAB-THESIS-2015-17


Expect verification with experimental data soon!




Real MicroBooNE PMT data





Real MicroBooNE PMT data





Real MicroBooNE PMT data




MicroBooNE system







TPB Efficiency



- WLS efficiency of pure TPB has been quantified – at 128nm, about 1.2 visible photons out for 1 UV photon in.




TPB coating efficiency



TPB time response

E. Segreto
Phys. Rev. C 91, 035503


- A 250 nm and above, TPB response is prompt (1-2 ns)
- Recent work shows that at 128 nm, TPB has a long time tail
- Comparison with scintillation profiles gives possible interpretation:
- Scintillating secondary electron liberated from the molecule, exciting TPB triplets
- This also helps explain the 1.2 shifting efficiency



	Decay time (ns)	Abundance (%)
Instantaneous component	1–10	60 ± 1
Intermediate component	49 ± 1	30 ± 1
Long component	3550 ± 500	8 ± 1
Spurious component	309 ± 10	2 ± 1

Environmental sensitivity of TPB

- We found that TPB is very sensitive to UV light and degrades in performance
- Also yellows in color with exposure
- Components in modern LArTPCs must last 10s of years



Journal of Instrumentation > Volume 7 > July 2012

C S Chiu et al 2012 JINST 7 P07007 doi:10.1088/1748-0221/7/07/P07007

Environmental effects on TPB wavelength-shifting coatings


C S Chiu, C Ignarra, L Bugel, H Chen, J M Conrad, B J P Jones, T Katori and I Moulton

Photodegradation Mechanism of TPB

NMR investigations show
spontaneously evolved photo-
initiating impurities in TPB


The impurity is present at low level
in lower commercially available
grades of TPB

Scintillation grade




Standard grade


*(decayed scint
grade shows
same peaks)*




1H selm1gp.2 7.04ppm selectively excited



1H selm1gp.2 7.73ppm selectively excited



Monitoring decay with GCMS



Working with GCMS we identified the degradation mechanism – radical mediated photo-oxidation to benzophenone

Demonstrated delayed degradation using radical mediators confirms this interpretation

Oxygen is required.
Submerged in LAr, should be safe from decay.

More on this later...

[Journal of Instrumentation > Volume 8 > January 2013](#)

B J P Jones *et al* 2013 *JINST* **8** P01013 doi:[10.1088/1748-0221/8/01/P01013](https://doi.org/10.1088/1748-0221/8/01/P01013)

Photodegradation mechanisms of tetraphenyl butadiene coatings for liquid argon detectors

B J P Jones^a, J K VanGemert^b, J M Conrad^a and A Pla-Dalmau^b

Bo Vertical Slice Test Stand


- Test bed for the MicroBooNE optical system including:
 - *Cryogenic photomultiplier tubes*
 - *Base electronics*
 - *Wavelength shifting plate*
 - *High voltage system + interlocks*
 - *Cables and splitters*
 - *Readout electronics*
 - *Cryostat feedthrough*
 - *Trace impurity monitors*
 - *Etc...*

Also used for liquid argon R&D work for MicroBooNE, DUNE, etc



Removing H₂O and O₂


- **Water** removed by 4A molecular sieves



*Very large effective surface area
for molecules with $d < 4\text{\AA}$*

water $d \sim 2.75\text{\AA}$

- **Oxygen** removed by activated copper granules




*Oxygen trapped by oxidizing
the copper*




- 10 ppt levels of water and oxygen impurity can be achieved

Effect of Nitrogen : 1. Quenching




*Competing Excimer
Dissociation Process*



*Rate dependent on the density of excimers and
density of impurity*

Quenching due to Nitrogen




Much stronger effect on triplet state – longer time for collisions to occur

Affects time constant as well as total light yield and fast / slow ratio

Systematic study shows that nitrogen should be kept < 2 ppm to prevent significant quenching

2010 JINST 5 P06003
WArP collaboration


Effect of Nitrogen: 2. Absorption




A measurement of the absorption of liquid argon scintillation light by dissolved nitrogen at the part-per-million level

B J P Jones, C S Chiu, J M Conrad, C M Ignarra, T Katori and M Touns

Published 24 July 2013 • © 2013 IOP Publishing Ltd and Sissa Medialab srl • [Journal of Instrumentation](#), Volume 8, July 2013



Top View






c1


File Edit View Options Tools

Help

Alpha Pulse Sizes



Implications for LArTPC Detectors




What should drive the design? (my opinion...)

MicroBooNE Event Zero:

(half design field)

Run 1147 Event 0. August 6th 2015 16:59

μ BooNE



What should drive the design? (my opinion...)

- **Maximize spatial resolution.**

What should drive the design? (my opinion...)

- **Maximize spatial resolution.**

Within the constraints:

- *Fixed \$\$*
- *Keep within detector real-estate*
- *Don't hurt the TPC system*

By doing the following:

- *More channels per \$*
- *More light per channel*

I am only going to highlight only the ideas that seem to achieve this.


More channels per \$

Denver Whittington –
LIDINE2015 talk



Light Guides for Large-Area Photon Detection



- Large active-area UV-collecting light guides
 - Acrylic or polystyrene imbued with wavelength-shifting compound
 - Based on design pioneered by MIT
 - Dip-coating w/ TPB in solvent (after studying many different methods)
 - 430 nm light propagated by total internal reflection to end




- Embed PD paddles inside anode plane assembly behind collection wires
 - Large photosensitive area with small photocathode area
 - Low-voltage SiPM bias
 - Easily scalable



- PMTs → SiPMs drives unit cost down
- Removes HV requirement
- But mandates additional collectors (eg lightguides)


More channels per \$:

This is an event
from the ICARUS
TPC:



More channels per \$:

This is an event
from the ICARUS
TPC:






Crosstalk from the PMT
system to the wire planes

Readout PMT system through crosstalk alone?

arxiv:1507.01997 Z. Moss et al

- Proposal to remove PMT readout electronics system altogether and purely read out via crosstalk
- Possibly boost crosstalk with an extra antenna to increase capacitance between systems
- No readout electronics means more \$ for sensitive elements




Back to light guides:

Denver Whittington –
LIDINE2015 talk



Light Guides for Large-Area Photon Detection



- Large active-area UV-collecting light guides
 - Acrylic or polystyrene imbued with wavelength-shifting compound
 - Based on design pioneered by MIT
 - Dip-coating w/ TPB in solvent (after studying many different methods)
 - 430 nm light propagated by total internal reflection to end



- Embed PD paddles inside anode plane assembly behind collection wires
 - Large photosensitive area with small photocathode area
 - Low-voltage SiPM bias
 - Easily scalable



- PMTs → SiPMs drives unit cost down
- Mandates additional collectors (eg lightguides)
- Collector requires optimization


Improved TPB Coating

Jarrett Moon – LIDINE2015 talk

- Increase TPB/acrylic ratio

0.5 g TPB	→	0.1 g TPB
50 mL toluene	→	50 mL toluene
10 mL ethanol	→	12 mL ethanol
1 g acrylic	→	0.1 g acrylic

- Initial test done in air, and showed *drastic* improvement





But for such a seemingly drastic improvement, the proof is in the pudding

(pudding = LAr)


Competing light guide designs

Denver Whittington – LIDINE2015 talk


Dip-Coated Acrylic (New MIT Recipe)




WLS Fibers + TPB Radiator



WLS Fibers in TPB-Coated Panel



TPB Dip-Coated Acrylic




Improvement by chemical stabilization

Improved initial performance and stabilization can be achieved with radical mediators

Interpretation: Partial TPB→BP decay in solution during coating preparation

More exploration of stabilizers could yield brighter TPB-PS coatings




[Journal of Instrumentation](#) > Volume 8 > January 2013

B J P Jones *et al* 2013 *JINST* **8** P01013 doi:[10.1088/1748-0221/8/01/P01013](https://doi.org/10.1088/1748-0221/8/01/P01013)


Photodegradation mechanisms of tetraphenyl butadiene coatings for liquid argon detectors

B J P Jones^a, J K VanGemert^b, J M Conrad^a and A Pla-Dalmau^b

LArIAT-style foils




Reflector foil before/after TPB evaporation



Test-mount of mock foil masks.


Will Foreman – LIDINE 2015 talk

Reflector-based solution (LArIAT)







Advantages of foils

- Detector is small, so cosmic background not problematic.
- LArIAT collection all light at one spot – foils homogenize detector response (good!)
- Detector goal is to study charge and light emission – so get lots of it (good!)
- On the flipside, very little spatial information in the reflected light.





The ARAPUCA

Ettore Segretto – LIDINE 2015 talk



Xe Doping

- Addition of xenon can produce excitation transfer which moves light to earlier times.
- **Common myth** : “TPB is brighter at LXe wavelengths” – not true!
- But bringing late light to earlier times after the pulse makes it more easily collectable, and more useful.
- ~100 ppm required → multiple tons for DUNE. Not a cheap proposal.



Summary

- Liquid argon TPCs are an important part of the US neutrino physics program
- Light collection systems are important for background rejection, drift-position measurement
- Significant R&D was required to realize present-generation experiments
- R&D for next-generation experiments should focus (IMO) on spatial specificity for geometrical tagging
- This can be achieved with more channels per \$ and more light per channel
- There are old ideas, new ideas, and ideas yet-to-be-had that can make this happen.